

RPCB2023

---

## Product Information

<b>Product SKU:</b>	RPCB2023	<b>Gene ID:</b>	13030	<b>Size:</b>	50µg
<b>Tag:</b>	C-His	<b>Reactivity:</b>	Mouse		

---

## Additional Information

<b>Expression Host:</b>	HEK293 cells	<b>Swissprot:</b>	P10605
<b>Purity:</b>	>95% as determined by SDS-PAGE		

---

## Protein Information

**Background:** -

**Protein Description:** High quality, high purity and low endotoxin recombinant Recombinant Mouse Cathepsin B Protein , tested reactivity in HEK293 cells and has been validated in SDS-PAGE.100% guaranteed.

**Endotoxin:** <1EU/µg

**Formulation:** Lyophilized from a 0.22 µm filtered solution of PBS, pH 7.4.

**Storage:** Store at -20°C.Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt.After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.