

Recombinant Human PPlase/FKBP7 Protein (aa 24-222, His Tag)

RPES1621

Description

This high-purity recombinant protein is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Protein Information

SKU: RPES1621

Purity: > 95 % as determined by reducing SDS-PAGE.

Contents: 50µg, 10µg
Bradford Reagent: 1 vial (2ml)

Concentration: Subject to label value.

Species: Human

Endotoxin: < 1.0 EU per µg of the protein as determined by the LAL method.

Synonyms: 23 kDa FK506-Binding Protein, 23 kDa FKBP, FK506-Binding Protein 7, FKBP-23, FKBP-7, FKBP23, FKBP7, PPlase FKBP7, Peptidyl-Prolyl Cis-Trans Isomerase FKBP7, Rotamase

Storage: Store at < -20°C, stable for 6 months. Please minimize freeze-thaw cycles. Store Bradford Reagent at Room Temperature for 1 year.

Tag: C-His

Shipping: This product is provided as liquid. It is shipped at frozen temperature with blue ice/gel packs. Upon receipt, store it immediately at < - 20°C.

Expression Host: HEK293 Cells

Bio-Activity: Not validated for activity

Calculated MW: 23.9 kDa

Formulation: Supplied as a 0.2 µm filtered solution of 20mM Tris-HCl, 150mM NaCl, 1mM CaCl₂ , 10% Glycerol, pH 7.5.

Observed MW: 25-32 kDa

Reconstitution: -

Accession: Q9Y680

Source: HEK293 Cells-derived Human PPlase, FKBP7 protein Gln24-Leu222, with an C-terminal His

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Manufacturers Statement: This final kit system is assembled and quality-released by Assay Genie Limited.

Sequence: Gln24-Leu222

Notes: Centrifuge before opening to ensure complete recovery of vial contents.

Form: Liquid