

Recombinant Mouse SerpinE2/SERPINE2 Protein (His Tag)

RPES1692

Description

This high-purity recombinant protein is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Protein Information

SKU: RPES1692

Purity: > 95 % as determined by reducing SDS-PAGE.

Contents: 50µg, 10µg
Bradford Reagent: 1 vial (2ml)

Concentration: Subject to label value.

Species: Mouse

Endotoxin: < 1.0 EU per µg of the protein as determined by the LAL method.

Synonyms: GDN, Glia-derived nexin, PI-7, PN-1, Peptidase inhibitor 7, Pi7, Pn1, Protease nexin 1, Protease nexin I, Serine protease-inhibitor 4, Serpin E2, Spi4

Storage: Store at < -20°C, stable for 6 months. Please minimize freeze-thaw cycles.
Store Bradford Reagent at Room Temperature for 1 year.

Tag: C-His

Shipping: This product is provided as liquid. It is shipped at frozen temperature with blue ice/gel packs. Upon receipt, store it immediately at < - 20°C.

Expression Host: HEK293 Cells

Bio-Activity: Not validated for activity

Calculated MW: 43.2 kDa

Formulation: Supplied as a 0.2 µm filtered solution of 20mM Sodium Acetate, 150 mM NaCl, 5% Trehalose, 5% Mannitol, 0.02% Tween 80, 1mM EDTA, pH 4.0.

Observed MW: 50-55 kDa

Reconstitution: -

Accession: Q07235

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Source: HEK293 Cells-derived Mouse SerpinE2/SERPINE2 protein Ser20-Pro397, with an C- terminal His

Manufacturers Statement: This final kit system is assembled and quality-released by Assay Genie Limited.

Sequence: Ser20-Pro397

Notes: Centrifuge before opening to ensure complete recovery of vial contents.

Form: Liquid