

# Recombinant Human PDIA5 Protein (His Tag)

RPES3106

## Description

---

This high-purity recombinant protein is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

## Protein Information

---

**SKU:** RPES3106

**Purity:** > 90 % as determined by reducing SDS-PAGE.

**Contents:** 50µg, 10µg  
Bradford Reagent: 1 vial (2ml)

**Concentration:** Subject to label value.

**Species:** Human

**Endotoxin:** < 1.0 EU per µg of the protein as determined by the LAL method.

**Synonyms:** PDIA5, PDIR, Protein disulfide isomerase-related protein, Protein disulfide-isomerase A5

**Storage:** Store at < -20°C, stable for 6 months. Please minimize freeze-thaw cycles.  
Store Bradford Reagent at Room Temperature for 1 year.

**Tag:** C-His

**Shipping:** This product is provided as liquid. It is shipped at frozen temperature with blue ice/gel packs. Upon receipt, store it immediately at < - 20°C.

**Expression Host:** HEK293 Cells

**Bio-Activity:** Not validated for activity

**Calculated MW:** 28.8 kDa

**Formulation:** Supplied as a 0.2 µm filtered solution of 20mM Acetate, 8% Sucrose, 0.05% Tween 80, pH5.0.

**Observed MW:** 30 kDa

**Reconstitution:** -

**Accession:** AAH01625.1

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

**Source:** HEK293 Cells-derived Human PDIA5 protein Ser22-Leu262, with an C-terminal His

**Sequence:** Ser22-Leu262

**Notes:** Centrifuge before opening to ensure complete recovery of vial contents.

**Manufacturers Statement:** This final kit system is assembled and quality-released by Assay Genie Limited.

**Form:** Liquid