

# Recombinant Human S100A6 Protein (E.coli, His Tag)

RPES3311

## Description

---

This high-purity recombinant protein is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

## Protein Information

---

**SKU:** RPES3311

**Purity:** > 95 % as determined by reducing SDS-PAGE.

**Contents:** 50µg, 10µg  
Bradford Reagent: 1 vial (2ml)

**Concentration:** Subject to label value.

**Species:** Human

**Endotoxin:** < 1.0 EU per µg of the protein as determined by the LAL method.

**Synonyms:** 2A9, 5B10, CABP, CACY, Calcyclin, Growth factor-inducible protein 2A9, MLN 4, PRA, Prolactin receptor-associated protein, Protein S100-A6, S100 calcium-binding protein A6, S100A6

**Storage:** Store at < -20°C, stable for 6 months. Please minimize freeze-thaw cycles.  
Store Bradford Reagent at Room Temperature for 1 year.

**Tag:** N-His

**Shipping:** This product is provided as liquid. It is shipped at frozen temperature with blue ice/gel packs. Upon receipt, store it immediately at < - 20°C.

**Expression Host:** E.coli

**Bio-Activity:** Not validated for activity

**Calculated MW:** 12.5 kDa

**Formulation:** Supplied as a 0.2 µm filtered solution of 20mM Tris-HCl, 150mM NaCl, 20% Glycerol, 1mM EDTA, pH 8.0.

**Observed MW:** 12 kDa

**Reconstitution:** -

**Accession:** P06703

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

**Source:** E.coli-derived Human S100A6 protein Met 1-Gly90, with an N-terminal His

**Manufacturers Statement:** This final kit system is assembled and quality-released by Assay Genie Limited.

**Sequence:** Met 1-Gly90

**Notes:** Centrifuge before opening to ensure complete recovery of vial contents.

**Form:** Liquid