

HDBS0160

Product Information

Product SKU:	HDBS0160	Host Species:	Homo sapiens	Size:	100 μ g
Target:	TNF α	Reactivity:	Human		

Additional Information

Conjugate:	Unconjugated	Uniprot:	P01375
Isotype:	IgG1		

Protein Information

Background: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Synonyms: DIF,TNF-alpha,TNFA,TNFSF2,cachexin,cachectin,TNF α

Common Name: CNTO 148

Recommended Dilutions: ELISA 1:5000-10000, Flow Cyt 1:100

Applications: ELISA, Flow Cyt

Usage: Research use only

Formulation: Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions.

Storage: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).Lyophilized antibodies are shipped at ambient temperature.