

Biotinylated Anti-TweakR (Enavatuzumab biosimilar) mAb

HDBS0247



Product Information

Product SKU: HDBS0247 **Host Species:** Humanized **Size:** 100 µg
Target: TweakR **Reactivity:** Human

Additional Information

Conjugate: Biotinylated **Uniprot:** Q9NP84
Isotype: IgG1

Protein Information

Background: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Synonyms: TNFRSF12A;FGF-inducible 14;FN14;TweakR;CD266

Common Name: PDL 192

Recommended Dilutions: ELISA 1:5000-10000; Flow Cyt 1:100

Applications: ELISA, Flow Cyt

Usage: Research use only

Formulation: Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

Storage: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.

Contact Details | Dublin, Ireland

Email: techsupport@assaygenie.com | Web: www.assaygenie.com

Copyright © 2024 Assay Genie Ltd, All Rights Reserved. All information / detail is correct at time of going to print.