

Biotinylated Anti-IFNG (Emapalumab biosimilar) mAb

HDBS0341



Product Information

Product SKU:	HDBS0341	Host Species:	Homo sapiens	Size:	100 µg
Target:	IFNG	Reactivity:	Human		

Additional Information

Conjugate:	Biotinylated	Uniprot:	P01579
Isotype:	IgG1		

Protein Information

Background: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Synonyms: IFN-gamma; Immune interferon

Common Name: NI-0501, emapalumab-lzsg

Recommended Dilutions: Flow Cyt 1:100

Applications: Flow Cyt

Usage: Research use only

Formulation: Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

Storage: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.

Contact Details | Dublin, Ireland

Email: techsupport@assaygenie.com | Web: www.assaygenie.com

Copyright © 2024 Assay Genie Ltd, All Rights Reserved. All information / detail is correct at time of going to print.