

HDBS0429

Product Information

Product SKU:	HDBS0429	Host Species:	Humanized	Size:	100 µg
Target:	CDH6	Reactivity:	Human		

Additional Information

Conjugate:	Unconjugated	Uniprot:	P55285
Isotype:	IgG1(E239D, M241L)		

Protein Information

Background: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Synonyms: CADH6;K-cadherin

Common Name: AMT-707; CUSP-06

Recommended Dilutions: ELISA 1:5000-10000, Flow Cyt 1:100

Applications: ELISA, Flow Cyt

Usage: Research use only

Formulation: Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

Storage: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized antibodies are shipped at ambient temperature.