

HDBS0722

---

**Product Information**

<b>Product SKU:</b>	HDBS0722	<b>Host Species:</b>	Whole mAb	<b>Size:</b>	100 µg
<b>Target:</b>	CTLA-4	<b>Reactivity:</b>	Human		

---

**Additional Information**

<b>Conjugate:</b>	Unconjugated	<b>Uniprot:</b>	P16410
<b>Isotype:</b>	IgG1		

---

**Protein Information**

**Background:** Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

**Synonyms:** CD;GSE;GRD4;ALPS5;CD152;CTLA-4;IDDM12;CELIAC3

**Common Name:** BMS-986249;BMS 986249;BMS986249

**Recommended Dilutions:** -

**Applications:** -

**Usage:** Research use only

**Formulation:** Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

**Storage:** Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.