

HDBS0791

Product Information

Product SKU: HDBS0791 **Host Species:** Whole mAb (Canine) **Size:** 100 µg
Target: Canine IL31 **Reactivity:** Canine

Additional Information

Conjugate: Unconjugated **Uniprot:** C7G0W1
Isotype: IgG2

Protein Information

Background: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Synonyms: Interleukin-31;IL-31

Common Name: ZTS-00008183;ZTS00008183;ZTS 00008183;ZTS-8183;ZTS8183;ZTS 8183

Recommended Dilutions: -

Applications: -

Usage: Research use only

Formulation: Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

Storage: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.