

# Biotinylated Anti-SIGLEC15 (Aldastotug biosimilar) mAb

HDBS0823



## Product Information

<b>Product SKU:</b>	HDBS0823	<b>Host Species:</b>	Whole mAb	<b>Size:</b>	100 µg
<b>Target:</b>	SIGLEC15	<b>Reactivity:</b>	Human		

## Additional Information

<b>Conjugate:</b>	Biotinylated	<b>Uniprot:</b>	Q6ZMC9
<b>Isotype:</b>	IgG1		

## Protein Information

**Background:** Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

**Synonyms:** CD33L3;HsT1361;SIGLEC-15

**Common Name:** BSI 060T;BSI-060T;BSI060T;PYX 106;PYX-106;PYX106

**Recommended Dilutions:** -

**Applications:** -

**Usage:** Research use only

**Formulation:** Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

**Storage:** Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.

Contact Details | Dublin, Ireland

Email: [techsupport@assaygenie.com](mailto:techsupport@assaygenie.com) | Web: [www.assaygenie.com](http://www.assaygenie.com)

Copyright © 2024 Assay Genie Ltd, All Rights Reserved. All information / detail is correct at time of going to print.