

HDBS1170

Product Information

Product SKU:	HDBS1170	Host Species:	Whole mAb	Size:	100 µg
Target:	INSR	Reactivity:	Human		

Additional Information

Conjugate:	Biotinylated	Uniprot:	P06213
Isotype:	IgG2		

Protein Information

**Background:** Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

**Synonyms:** HHF5;CD220

**Common Name:** RZ-358;X 358;XOMA 247;XOMA-247;XOMA247;XOMA 358;XOMA-358;XOMA358;XPA.15.247

**Recommended Dilutions:** -

**Applications:** -

**Usage:** Research use only

**Formulation:** Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

**Storage:** Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.