

# Anti-Human CD257 (BAFF) (Tabalumab) - Biotin Biosimilar

IVMB0441

## Description

---

This Anti-Human CD257 (BAFF) (Tabalumab) - Biotin Biosimilar is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

## Product Information

---

**SKU:** IVMB0441

**Contents:** 100ug  
Bradford Reagent: 1 vial (2ml)

**Synonyms:** Tabalumab, CD257, BAFF, TNFSF13b, BLYS

**Category:** Biosimilar Recombinant Human Monoclonal Antibody

**Target:** BAFF

**Clone:** Tabalumab

**Isotype:** Human IgG1κ

**Applications:** -

**Specificity:** This non-therapeutic biosimilar antibody uses the same variable region sequence as the therapeutic antibody Tabalumab. Tabalumab neutralizes soluble human, cynomolgus monkey, and rabbit BAFF. Additionally, Tabalumab neutralizes membrane-bound BAFF. This product is for research use only.

## Antibody Data

---

**Reactivity:** Cynomolgus Monkey, Human, Rabbit

**Host species:** Human

**Expression Host:** HEK-293 Cells

### Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

<b>Immunogen:</b>	Original Antibody was raised against soluble Human BAFF.
<b>Product concentration:</b>	0.5 mg/ml
<b>Endotoxin Level:</b>	-
<b>Purity:</b>	-
<b>Formulation:</b>	This Biotinylated antibody is formulated in 0.01 M phosphate buffered saline (150 mM NaCl) PBS pH 7.4, 1% BSA and 0.09% sodium azide as a preservative.

## Preparation & Storage

---

<b>Storage:</b>	This biotinylated antibody is stable when stored at 2-8°C. Do not freeze. Store Bradford Reagent at Room Temperature for 1 Year.
<b>Shipping:</b>	2-8°C Wet Ice
<b>Preparation:</b>	-
<b>Recommended Dilution Buffer:</b>	-
<b>Recommended Usage:</b>	-
<b>Protein Quantification (Optional):</b>	To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <a href="https://www.assaygenie.com/bradford-protein-assay-protocol/">https://www.assaygenie.com/bradford-protein-assay-protocol/</a> to view the full protocol