

# VivoGenie Anti-Human RANKL (Denosumab) Biosimilar

IVMB0490

## Description

---

This VivoGenie Anti-Human RANKL (Denosumab) Biosimilar is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

## Product Information

---

<b>SKU:</b>	IVMB0490
<b>Contents:</b>	1mg, 5mg, 25mg, 50mg, 100mg Bradford Reagent: 1 vial (2ml)
<b>Synonyms:</b>	Osteoprotegerin Ligand (OPGL), Osteoclast Differentiation Factor (ODF), TNF related Activation-Induced Cytokine (TRANCE), Tumor Necrosis Factor Ligand Superfamily Member 11 (TNFSF11)
<b>Category:</b>	Biosimilar Recombinant Human Monoclonal Antibody
<b>Target:</b>	RANKL
<b>Clone:</b>	AMG162
<b>Isotype:</b>	Human IgG2k
<b>Applications:</b>	<a href="#">ELISA</a> <a href="#">FA</a> <a href="#">IHC</a> <a href="#">WB</a>
<b>Specificity:</b>	This non-therapeutic biosimilar antibody uses the same variable region sequence as the therapeutic antibody Denosumab. This product is for research use only. Denosumab activity is directed against human RANKL (receptor activator of NFkB ligand).

## Antibody Data

---

<b>Reactivity:</b>	Human
<b>Host species:</b>	Human

### Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

<b>Expression Host:</b>	HEK-293 Cells
<b>Immunogen:</b>	Purified Recombinant Human RANKL
<b>Product concentration:</b>	≥ 5.0 mg/ml
<b>Endotoxin Level:</b>	< 1.0 EU/mg as determined by the LAL method
<b>Purity:</b>	≥95% by SDS Page, ≥95% Monomer by analytical SEC
<b>Formulation:</b>	This biosimilar antibody is aseptically packaged and formulated in 0.01 M phosphate buffered saline (150 mM NaCl) PBS pH 7.2 - 7.4 with no carrier protein, potassium, calcium or preservatives added. Due to inherent biochemical properties of antibodies, certain products may be prone to precipitation over time. Precipitation may be removed by aseptic centrifugation and/or filtration.

## Preparation & Storage

---

<b>Storage:</b>	Functional grade preclinical antibodies may be stored sterile as received at 2-8°C for up to one month. For longer term storage, aseptically aliquot in working volumes without diluting and store at ≤ -70°C. Avoid Repeated Freeze Thaw Cycles. Store Bradford Reagent at Room Temperature for 1 Year.
<b>Shipping:</b>	2-8°C Wet Ice
<b>Preparation:</b>	Recombinant biosimilar antibodies are manufactured in an animal free facility using only in vitro protein free cell culture techniques and are purified by a multi-step process including the use of protein A or G to assure extremely low levels of endotoxins, leachable protein A or aggregates.
<b>Recommended Dilution Buffer:</b>	-
<b>Recommended Usage:</b>	-
<b>Protein Quantification (Optional):</b>	To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <a href="https://www.assaygenie.com/bradford-protein-assay-protocol/">https://www.assaygenie.com/bradford-protein-assay-protocol/</a> to view the full protocol