



Product Datasheet

Purified Anti-Mouse MHC II (I-A/I-E) Antibody [M5/114] Catalogue Code: AGEL0018

Antibody Data

Product SKU:	AGEL0018	Clone:	M5/114
Applications:	FCM		
Reactivity:	Mouse		

Important Note:

Centrifuge before opening to ensure complete recovery of vial contents.

Product Information:

Alternate Names:	H2-Ab1/Eb1;Major histocompatibility protein class II beta chain;MHC class II H2-IA-beta-psi;l-E beta MHC class II;MHC class II;
Uniprot ID:	P14483 O78196
Background:	These class II molecules are expressed on antigen presenting cells (including B cells) and a subset of T cells from H-2b,d,q,r bearing mice and are involved in antigen presentation to T cells expressing CD3/TCR and CD4 proteins.
Form:	Liquid
Conjugation:	Unconjugated
Size:	25µg, 100µg
Host Species:	Rat
Isotype:	Rat IgG2b, κ
Isotype Control:	Purified Rat IgG2b, κ Isotype Control[LTF-2] [Product AGEL0018]
Storage Buffer:	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer.
Shipping:	Biological ice pack at 4°C

Stability & Storage: Keep as concentrated solution. Store at 2~8°C and protected from prolonged exposure to light. Do not freeze. Centrifuge before opening to ensure complete recovery of vial contents. This product is guaranteed up to one year from purchase.

Recommended Usage: Each lot of this antibody is quality control tested by flow cytometric analysis. For flow cytometric staining, the suggested use of this reagent is $\leq 0.25 \mu\text{g}$ per 10^6 cells in 100 μL volume or 100 μL of whole blood. It is recommended that the reagent be titrated for optimal performance for each application.