



## Product Datasheet

### Purified Anti-Mouse CD48 Antibody [HM48-1]

Catalogue Code: AGEL0109

#### Antibody Data

<b>Product SKU:</b>	<b>AGEL0109</b>	<b>Clone:</b>	<b>HM48-1</b>
<b>Applications:</b>	<b>FCM</b>		
<b>Reactivity:</b>	<b>Mouse</b>		

#### Important Note:

Centrifuge before opening to ensure complete recovery of vial contents.

#### Product Information:

<b>Alternate Names:</b>	CD48 antigen;Cd48;BCM1 surface antigen;BLAST-1;HM48-1;MRC OX-45 surface antigen;SLAMF2;sgp-60;CD48;
<b>Uniprot ID:</b>	P18181
<b>Background:</b>	CD48 is a 45 kD GPI-anchored glycoprotein also known as BCM1, Blast-1 (human), and OX-45 (rat). It is a member of the Ig superfamily, expressed on T and B cells and monocytes/macrophages. It plays a role in adhesion and T cell recognition. The primary ligands for CD48 are CD2 and CD244.
<b>Form:</b>	Liquid
<b>Conjugation:</b>	Unconjugated
<b>Size:</b>	25&micro;g, 100&micro;g
<b>Host Species:</b>	Armenian Hamster
<b>Isotype:</b>	Armenian Hamster IgG
<b>Isotype Control:</b>	Purified Armenian Hamster IgG Isotype Control[PIP] [Product AGEL0109]
<b>Storage Buffer:</b>	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer.
<b>Shipping:</b>	Biological ice pack at 4°C

**Stability & Storage:** Keep as concentrated solution. Store at 2~8°C and protected from prolonged exposure to light. Do not freeze. Centrifuge before opening to ensure complete recovery of vial contents. This product is guaranteed up to one year from purchase.

**Recommended Usage:** Each lot of this antibody is quality control tested by flow cytometric analysis. For flow cytometric staining, the suggested use of this reagent is  $\leq 0.25 \mu\text{g}$  per  $10^6$  cells in 100  $\mu\text{L}$  volume or 100  $\mu\text{L}$  of whole blood. It is recommended that the reagent be titrated for optimal performance for each application.

---