

**Low Endotoxin Purified Anti-Mouse  
F4/80 Antibody [Cl:A3-1]**

Catalogue Code: AGEL0561

**Antibody Data**

<b>Product SKU:</b>	<b>AGEL0561</b>	<b>Clone:</b>	<b>Cl:A3-1</b>
<b>Applications:</b>	<b>FCM;Depletion;FA</b>		
<b>Reactivity:</b>	<b>Mouse</b>		

**Important Note:**

Centrifuge before opening to ensure complete recovery of vial contents.

**Product Information:**

<b>Alternate Names:</b>	Adhesion G protein-coupled receptor E1;Adgre1;Cell surface glycoprotein F4/80;EGF-like module receptor 1;Adgre1;Emr1; Gpf480;
<b>Uniprot ID:</b>	Q61549
<b>Background:</b>	F4/80 is a 160 kD glycoprotein. It is characterized as a member of the epidermal growth factor (EGF)-transmembrane 7 (TM7) family. F4/80, also known as EMR1 or Ly71, has been widely used as a murine macrophage marker, which is expressed on the majority of tissue macrophages including peritoneal macrophages, macrophages in lung, gut, thymus and red pulp of spleen (but not on the macrophages located in T cell areas of the spleen, lymph node and Peyer's patch), Kuffer cells, Langerhans cells, and bone marrow stromal cells. F4/80 has also been shown on a subset of dendritic cells. The biological ligand of F4/80 has not been identified, but it has been reported that F4/80 is required for induction of CD8+ T cells-mediated peripheral tolerance.
<b>Form:</b>	Liquid
<b>Conjugation:</b>	None (AF/LE)
<b>Size:</b>	50&micro;g, 500&micro;g, 1mg
<b>Host Species:</b>	Rat
<b>Isotype:</b>	Rat IgG2b, $\kappa$
<b>Isotype Control:</b>	AF/LE Purified Rat IgG2b, $\kappa$ Isotype Control[LTF-2] [Product AGEL0561]
<b>Storage Buffer:</b>	0.2 $\mu$ m filtered in PBS, pH 7.2. Azide Free (AF)/Low Endotoxin (LE): Contains no stabilizers or stabilizers. Endotoxin level is < 2 EU/mg as Determined by LAL gel clotting assay.
<b>Shipping:</b>	Biological ice pack at 4°C

**Stability & Storage:** Keep as concentrated solution. Store at 2~8°C and protected from prolonged exposure to light. Do not freeze. Centrifuge before opening to ensure complete recovery of vial contents. This product is guaranteed up to one year from purchase.

**Recommended Usage:** Each lot of this antibody is quality control tested by flow cytometric analysis. For flow cytometric staining, the suggested use of this reagent is  $\leq 1.0 \mu\text{g}$  per  $10^6$  cells in 100  $\mu\text{L}$  volume or 100  $\mu\text{L}$  of whole blood. It is recommended that the reagent be titrated for optimal performance for each application.