

Antibody Data

Product SKU:	AGEL1143	Clone:	TY/23
Applications:	FCM		
Reactivity:	Mouse		


Important Note:

Centrifuge before opening to ensure complete recovery of vial contents.

Product Information:

Alternate Names: 5'-nucleotidase;Nt5e;5'-NT;Ecto-5'-nucleotidase;CD73;Nt5; Nte;Nt5e;
Uniprot ID: Q61503
Background: CD73 (ecto-5`-nucleotidase) is a 69 kD GPI-anchored surface protein. In mice, expression of CD73 in bone marrow is restricted to CD11b+ myeloid cells. In spleen, it is largely expressed on T cells.

Form: Liquid
Conjugation: PE
Size: 25µg, 100µg
Host Species: Rat
Isotype: Rat IgG2a, κ



Isotype Control: PE Rat IgG2a, κ Isotype Control[2A3] [Product AGEL1143]
Storage Buffer: Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Shipping: Biological ice pack at 4°C

Stability & Storage: Keep as concentrated solution. Store at 2~8°C and protected from prolonged exposure to light. Do not freeze. Centrifuge before opening to ensure complete recovery of vial contents. This product is guaranteed up to one year from purchase.

Recommended Usage: Each lot of this antibody is quality control tested by flow cytometric analysis. Please check your vial before the experiment. Since applications vary, the appropriate dilutions must be determined for individual use. We suggest each investigator should titrate the reagent to obtain optimal results [The recommended concentration is 0.1-1 µg/10⁶ cells in 100 µL volume].
