



## Product Datasheet

# Low Endotoxin Purified Anti-Human CD44 Antibody [Hermes-1]

Catalogue Code: AGEL2101

### Antibody Data

<b>Product SKU:</b>	<b>AGEL2101</b>	<b>Clone:</b>	<b>Hermes-1</b>
<b>Applications:</b>	<b>FCM;Block</b>		
<b>Reactivity:</b>	<b>Human</b>		

### Important Note:

Centrifuge before opening to ensure complete recovery of vial contents.

### Product Information:

<b>Alternate Names:</b>	LHR;PGP-I;CDw44;Epican;HUTCH-I;
<b>Uniprot ID:</b>	P16070
<b>Background:</b>	CD44 is a 80-95 kD glycoprotein also known as Hermes, Pgp1, H-CAM, or HUTCH. It is expressed on all leukocytes, endothelial cells, hepatocytes, and mesenchymal cells. As B and T cells become activated or progress to the memory stage, CD44 expression increases from a low or mid level of intensity to high expression levels. Thus, CD44 has been reported to be a valuable marker for memory cell subsets. CD44 is an adhesion molecule involved in leukocyte attachment to and rolling on endothelial cells, homing to peripheral lymphoid organs and to the sites of inflammation, and leukocyte aggregation.
<b>Form:</b>	Liquid
<b>Conjugation:</b>	None (AF/LE)
<b>Size:</b>	50µg, 500µg, 1mg
<b>Host Species:</b>	Rat
<b>Isotype:</b>	Rat IgG2a, κ
<b>Isotype Control:</b>	AF/LE Purified Rat IgG2a, κ Isotype Control[2A3] [Product AGEL2101]
<b>Storage Buffer:</b>	0.2 µm filtered in PBS, pH 7.2. Azide Free (AF)/Low Endotoxin (LE): Contains no stabilizers or stabilizers. Endotoxin level is < 2 EU/mg as Determined by LAL gel clotting assay.
<b>Shipping:</b>	Biological ice pack at 4°C

**Stability & Storage:** Keep as concentrated solution. Store at 2~8°C and protected from prolonged exposure to light. Do not freeze. Centrifuge before opening to ensure complete recovery of vial contents. This product is guaranteed up to one year from purchase.

**Recommended Usage:** Each lot of this antibody is quality control tested by flow cytometric analysis. For flow cytometric staining, the suggested use of this reagent is  $\leq 2.0 \mu\text{g}$  per  $10^6$  cells in 100  $\mu\text{L}$  volume or 100  $\mu\text{L}$  of whole blood. It is recommended that the reagent be titrated for optimal performance for each application.