

**Low Endotoxin Purified Anti-Mouse
CD272 Antibody [PK18.6]**
Catalogue Code: AGEL2159

Antibody Data

| | | | |
|----------------------|-----------------|---------------|---------------|
| Product SKU: | AGEL2159 | Clone: | PK18.6 |
| Applications: | FCM | | |
| Reactivity: | Mouse | | |

Important Note:

Centrifuge before opening to ensure complete recovery of vial contents.

Product Information:

| | |
|-------------------------|---|
| Alternate Names: | B- and T-lymphocyte attenuator;Btla;B- and T-lymphocyte-associated protein;CD272;Btla; |
| Uniprot ID: | Q7TSA3 |
| Background: | CD272, also known as B and T lymphocyte attenuator (BTLA), is an Ig superfamily co-inhibitory receptor with structural similarity to programmed cell death 1 (PD-1) and CTLA-4. BTLA is expressed on B cells, T cells, macrophages, dendritic cells, NKT cells, and NK cells. Engagement of BTLA by its ligand herpes virus entry mediator (HVEM) is critical for negatively regulating immune response. The absence of BTLA with HVEM inhibitory interactions leads to increased experimental autoimmune encephalomyelitis severity, enhanced rejection of partially mismatched allografts, an increased CD8+ memory T cell population, increased severity of colitis, and reduced effectiveness of T regulatory cells. BTLA plays an important role in the induction of peripheral tolerance of both CD4+ and CD8+ T cells in vivo. Tolerant T cells have significantly higher expression of BTLA compared with effectors and naïve T cells. BTLA may cooperate with CTLA-4 and PD-1 to control T cell tolerance and autoimmunity. It was reported that BTLA may regulate T cell function by binding to B7-H4, but further studies are needed to confirm. The existence of three distinct BTLA alleles has been reported. |
| Form: | Liquid |
| Conjugation: | None (AF/LE) |
| Size: | 50µg, 500µg, 1mg |
| Host Species: | Rat |
| Isotype: | Rat IgG1, κ |
| Isotype Control: | AF/LE Purified Rat IgG1, κ Isotype Control[HRPN] [Product AGEL2159] |
| Storage Buffer: | 0.2 µm filtered in PBS, pH 7.2. Azide Free (AF)/Low Endotoxin (LE): Contains no stabilizers or stabilizers. Endotoxin level is < 2 EU/mg as Determined by LAL gel clotting assay. |
| Shipping: | Biological ice pack at 4°C |

Stability & Storage: Keep as concentrated solution. Store at 2~8°C and protected from prolonged exposure to light. Do not freeze. Centrifuge before opening to ensure complete recovery of vial contents. This product is guaranteed up to one year from purchase.

Recommended Usage: Each lot of this antibody is quality control tested by flow cytometric analysis. For flow cytometric staining, the suggested use of this reagent is $\leq 0.5 \mu\text{g}$ per 10^6 cells in 100 μL volume or 100 μL of whole blood. It is recommended that the reagent be titrated for optimal performance for each application.
