

Antibody Data

Product SKU:	AGEL2852	Clone:	PIP
Applications:	FCM		
Reactivity:	-		


Important Note:

Centrifuge before opening to ensure complete recovery of vial contents.

Product Information:

Alternate Names: -
Uniprot ID: -
Background: -

Form: Liquid
Conjugation: PE/Genie Fluor594
Size: 25µg, 100µg
Host Species: Armenian Hamster
Isotype: Armenian Hamster IgG



Isotype Control: -
Storage Buffer: Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Shipping: Biological ice pack at 4°C

Stability & Storage: Keep as concentrated solution. Store at 2~8°C and protected from prolonged exposure to light. Do not freeze. Centrifuge before opening to ensure complete recovery of vial contents. This product is guaranteed up to one year from purchase.

Recommended Usage: Each lot of this antibody is quality control tested by flow cytometric analysis as negative control. Use at concentrations comparable to those of the specific antibody of interest.