



Product Datasheet
Mouse IgG1 Isotype control monoclonal antibody
Catalogue Code: AGIM0104

Antibody Data

Product SKU:	AGIM0104	Clone	PPV-06
Applications:	FC	Conjugate	Unconjugated
Reactivity:	Mouse		

Immunogen Information:

Host Species: Mouse

Isotype: IgG1

Storage: Aliquot and store at -20°C. Avoid freeze/thaw cycles. Should be handled under aseptic conditions.

Background:

The specificity of staining by monoclonal antibodies to target antigens should be verified by establishing the amount of non-specific antibody binding. In general, non-reactive immunoglobulin of the same iso type is included as a negative control for each specific monoclonal antibody used in a particular immunoassay.