



Product Datasheet
Mouse IgG1 Isotype control monoclonal antibody (Biotin)
Catalogue Code: AGIM0105

Antibody Data

Product SKU:	AGIM0105	Clone	PPV-06
Applications:	FC	Conjugate	Biotin
Reactivity:	Mouse		

Immunogen Information:

Host Species: Mouse

Isotype: IgG1

Storage: Aliquot and store at -20°C. Avoid freeze/thaw cycles. Should be handled under aseptic conditions. Avoid freeze-thaw cycles

Background:

The specificity of staining by monoclonal antibodies to target antigens should be verified by establishing the amount of non-specific antibody binding. In general, non-reactive immunoglobulin of the same iso type is included as a negative control for each specific monoclonal antibody used in a particular immunoassay.