



Product Datasheet
Mouse IgG2a Isotype control monoclonal antibody (Biotin)
Catalogue Code: AGIM0110

Antibody Data

Product SKU:	AGIM0110	Clone	PPV-04
Applications:	FC	Conjugate	Biotin
Reactivity:	Mouse		

Immunogen Information:

Host Species: Mouse

Isotype: IgG2a

Storage: Store at 4 °C. For long term storage aliquot and store at 20°C. Avoid freeze/thaw cycles.

Background:

The specificity of staining by monoclonal antibodies to target antigens should be verified by establishing the amount of non-specific antibody binding. In general, non-reactive immunoglobulin of the same isotype is included as a negative control for each specific monoclonal antibody used in a particular immunoassay.