



Product Datasheet
Mouse IgG2b Isotype control monoclonal antibody
Catalogue Code: AGIM0114

Antibody Data

Product SKU:	AGIM0114	Clone	PLRV219
Applications:	FC	Conjugate	Unconjugated
Reactivity:	Mouse		

Immunogen Information:

Host Species: Mouse

Isotype: IgG2b

Storage: Store at 4°C. For long-term, aliquot and store at 20°C. Avoid freeze/thaw cycles.

Background:

The specificity of staining by monoclonal antibodies to target antigens should be verified by establishing the amount of non-specific antibody binding. In general, non-reactive immunoglobulin of the same isotype is included as a negative control for each specific monoclonal antibody used in a particular immunoassay.