



Product Datasheet  
**Mouse IgG2b control monoclonal antibody  
(FITC)**  
Catalogue Code: AGIM0116

## Antibody Data

<b>Product SKU:</b>	AGIM0116	<b>Clone</b>	PLRV219
<b>Applications:</b>	FC	<b>Conjugate</b>	FITC
<b>Reactivity:</b>	Mouse		

## Immunogen Information:

**Host Species:** Mouse

**Isotype:** IgG2b

**Storage:** Store at 4°C. Do not freeze. Avoid prolonged exposure to light.

## Background:

The specificity of staining by monoclonal antibodies to target antigens should be verified by establishing the amount of non-specific antibody binding. In general, non-reactive immunoglobulin of the same isotype is included as a negative control for each specific monoclonal antibody used in a particular immunoassay.