



Product Datasheet
**Mouse IgG2b control monoclonal antibody
(PE)**
Catalogue Code: AGIM0138

Antibody Data

| | | | |
|----------------------|----------|------------------|-------|
| Product SKU: | AGIM0138 | Clone | GC198 |
| Applications: | FC | Conjugate | PE |
| Reactivity: | Mouse | | |

Immunogen Information:

Host Species: Mouse

Isotype: IgG2b

Storage: Store at 4 °C. Do not freeze. Avoid prolonged exposure to light.

Background:

The specificity of staining by monoclonal antibodies to target antigens should be verified by establishing the amount of non-specific antibody binding. In general, non-reactive immunoglobulin of the same isotype is included as a negative control for each specific monoclonal antibody used in a particular immunoassay.