



Product Datasheet  
**Rat CD11a monoclonal antibody (FITC)**  
Catalogue Code: AGIM0279

### Antibody Data

<b>Product SKU:</b>	AGIM0279	<b>Clone</b>	WT.1
<b>Applications:</b>	FC	<b>Conjugate</b>	FITC
<b>Reactivity:</b>	Rat		

### Immunogen Information:

**Host Species:** Mouse

**Isotype:** IgG2a

**Storage:** Store at 4 °C. Avoid prolonged exposure to light.

### Background: