



Product Datasheet  
**Rat CD11b/c monoclonal antibody**  
Catalogue Code: AGIM0297

## Antibody Data

|                      |          |                  |              |
|----------------------|----------|------------------|--------------|
| <b>Product SKU:</b>  | AGIM0297 | <b>Clone</b>     | MRC OX-42    |
| <b>Applications:</b> | FC, IHC  | <b>Conjugate</b> | Unconjugated |
| <b>Reactivity:</b>   | Rat      |                  |              |

## Immunogen Information:

**Host Species:** Mouse

**Isotype:** IgG2a

**Storage:** For long term storage, aliquot and freeze unused portion at 20°C in volumes appropriate for single usage. Avoid freeze/thaw cycles.

## Background: