



Product Datasheet
Mouse IgG2b Isotype control monoclonal antibody (FITC)
Catalogue Code: AGIM0337

Antibody Data

Product SKU:	AGIM0337	Clone	TG1.7
Applications:	FC	Conjugate	FITC
Reactivity:	Mouse		

Immunogen Information:

Host Species: Mouse

Isotype: IgG2b

Storage: Store at 4 °C. Avoid prolonged exposure to light.

Background:

The specificity of staining by monoclonal antibodies to target antigens should be verified by establishing the amount of non-specific antibody binding. In general, non-reactive immunoglobulin of the same isotype is included as a negative control for each specific monoclonal antibody used in a particular immunoassay.