



Product Datasheet
Mouse IgG2a Isotype control monoclonal antibody (PerCP)
Catalogue Code: AGIM0396

Antibody Data

Product SKU:	AGIM0396	Clone	X5563
Applications:	FC	Conjugate	PerCP
Reactivity:	Mouse		

Immunogen Information:

Host Species: Mouse

Isotype: IgG2a

Storage: Store at 4 °C. Do not freeze. Avoid prolonged exposure to light.

Background:

The specificity of staining by monoclonal antibodies to target antigens should be verified by establishing the amount of non-specific antibody binding. In general, non-reactive immunoglobulin of the same isotype is included as a negative control for each specific monoclonal antibody used in a particular immunoassay.