

## CAB0536

---

### Product Information

<b>Product SKU:</b>	CAB0536	<b>Gene ID:</b>	2963	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	ARC2513	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human,Mouse,Rat

---

### Additional Information

<b>Observed MW:</b>	30kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	28kDa	<b>Isotype:</b>	IgG

---

### Immunogen Information

**Background:** Predicted to enable RNA polymerase II general transcription initiation factor activity. Involved in transcription by RNA polymerase II. Located in microtubule cytoskeleton and nucleoplasm. Part of transcription preinitiation complex.

**Recommended Dilution:** WB,1:500 - 1:1000 IHC-P,1:50 - 1:200

**Synonyms:** BTF4; RAP30; TF2F2; TFIIIF; GTF2F2

**Purification Method:** Affinity purification

**Immunogen:** A synthetic peptide corresponding to a sequence within amino acids 1-100 of human GTF2F2 (P13984).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.