

Product Information

Product SKU:	CAB0647	Gene ID:	6627	Size:	20uL, 100uL
Clone No:	ARC2532	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	28kDa	Conjugate:	Unconjugated
Calculated MW:	28kDa	Isotype:	IgG

Immunogen Information

**Background:** Enables RNA binding activity. Involved in mRNA splicing, via spliceosome and spermatogenesis. Located in nuclear speck. Part of U2-type catalytic step 2 spliceosome and U2-type precatalytic spliceosome. Implicated in connective tissue disease.

Recommended Dilution:	WB,1:500 - 1:1000 IHC-P,1:50 - 1:200
Synonyms:	Lea1; U2A'; SNRPA1
Purification Method:	Affinity purification
Immunogen:	A synthetic peptide corresponding to a sequence within amino acids 156-255 of human SNRPA1 (P09661).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.