

CAB0730

Product Information

Product SKU:	CAB0730	Gene ID:	7301	Size:	20uL, 100uL
Clone No:	ARC1841	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	130kDa	Conjugate:	Unconjugated
Calculated MW:	97kDa	Isotype:	IgG

Immunogen Information

Background: The gene is part of a 3-member transmembrane receptor kinase receptor family with a processed pseudogene distal on chromosome 15. The encoded protein is activated by the products of the growth arrest-specific gene 6 and protein S genes and is involved in controlling cell survival and proliferation, spermatogenesis, immunoregulation and phagocytosis. The encoded protein has also been identified as a cell entry factor for Ebola and Marburg viruses.

Recommended Dilution: WB,1:500 - 1:1000

Synonyms: BYK; Dtk; RSE; Rek; Sky; Tif; Etk-2; Tyro3

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 791-890 of human Tyro3 (Q06418).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.