

CAB0865

Product Information

Product SKU:	CAB0865	Gene ID:	56963	Size:	20uL, 100uL
Clone No:	ARC2546	Host Species:	Rabbit	Reactivity:	Mouse,Rat

Additional Information

Observed MW:	49kDa	Conjugate:	Unmodified
Calculated MW:	49kDa	Isotype:	IgG

Immunogen Information

Background: This gene encodes a member of the repulsive guidance molecule family. The encoded protein is a glycosylphosphatidylinositol-anchored glycoprotein that functions as an axon guidance protein in the developing and adult central nervous system. This protein may also function as a tumor suppressor in some cancers. Alternate splicing results in multiple transcript variants.

Recommended Dilution: WB,1:500 - 1:1000

Synonyms: RGM; RGMA

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 350-450 of human RGMA (Q96B86).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.

Contact Details | Dublin, Ireland

Email: techsupport@assaygenie.com | **Web:** www.assaygenie.com

Copyright © 2020 Reagent Genie, All Rights Reserved. All information / detail is correct at time of going to print.