

**CAB0982**

## Product Information

<b>Product SKU:</b>	CAB0982	<b>Gene ID:</b>	5291	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	-	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human,Mouse,Rat

## Additional Information

<b>Observed MW:</b>	110kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	123kDa	<b>Isotype:</b>	IgG

## Immunogen Information

**Background:** This gene encodes an isoform of the catalytic subunit of phosphoinositide 3-kinase (PI3K). These kinases are important in signaling pathways involving receptors on the outer membrane of eukaryotic cells and are named for their catalytic subunit. The encoded protein is the catalytic subunit for PI3Kbeta (PI3KB). PI3KB has been shown to be part of the activation pathway in neutrophils which have bound immune complexes at sites of injury or infection. Alternative splicing results in multiple transcript variants.

<b>Recommended Dilution:</b>	WB,1:1000 - 1:5000
<b>Synonyms:</b>	PI3K; PIK3C1; P110BETA; PI3KBETA; PI3 Kinase p110 beta
<b>Purification Method:</b>	Affinity purification
<b>Immunogen:</b>	Recombinant fusion protein containing a sequence corresponding to amino acids 410-600 of human PI3 Kinase p110 beta (NP_006210.1).
<b>Storage:</b>	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.