

EV71 3C Rabbit Polyclonal Antibody



CAB10003

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to Figures

Calculated MW:

Applications:

WB

Reactivity:

Enterovirus

Antibody Information

Recommended dilutions:

WB 1:1000 - 1:5000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

Uniprot

Synonyms:

Immunogen:

Recombinant protein of Enterovirus EV71 3C

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
