

CAB10003

Product Information

Product SKU:	CAB10003	Gene ID:	-	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Enterovirus

Additional Information

Observed MW:	Refer to Figures	Conjugate:	Unconjugated
Calculated MW:	-	Isotype:	IgG

Immunogen Information

Background:

Recommended Dilution:	WB, 1:1000 - 1:5000
Synonyms:	-
Purification Method:	Affinity purification
Immunogen:	Recombinant protein of Enterovirus EV71 3C
Storage:	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.