

**CAB10036**

---

**Product Information**

<b>Product SKU:</b>	CAB10036	<b>Gene ID:</b>	63893	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	-	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human,Mouse,Rat

---

**Additional Information**

<b>Observed MW:</b>	200kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	141kDa	<b>Isotype:</b>	IgG

---

**Immunogen Information**

**Background:** Enables ubiquitin conjugating enzyme activity and ubiquitin protein ligase activity. Involved in positive regulation of BMP signaling pathway; protein ubiquitination; and retrograde transport, endosome to Golgi. Located in cytoplasm and nuclear body.

**Recommended Dilution:** WB,1:1000 - 1:4000 IF/ICC,1:50 - 1:200 IP,0.5µg-4µg antibody for 200µg-400µg extracts of whole cells

**Synonyms:** E2-230K; 2O

**Purification Method:** Affinity purification

**Immunogen:** Recombinant fusion protein containing a sequence corresponding to amino acids 1073-1292 of human UBE2O (NP\_071349.3).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.