

LY75 Rabbit Polyclonal Antibody

CAB10065



Product Information

Product SKU:	CAB10065	Gene ID:	4065	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse

Additional Information

Observed MW:	250kDa	Conjugate:	Unconjugated
Calculated MW:	198kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to enable signaling receptor activity. Predicted to be involved in immune response and inflammatory response. Located in extracellular exosome.

Recommended Dilution:	WB,1:500 - 1:1000
Synonyms:	CD205; LY-75; CLEC13B; DEC-205; GP200-MR6; LY75
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 216-501 of human LY75 (NP_002340.2).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.