

LY75 Rabbit Polyclonal Antibody



CAB10065

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

205kDa

Calculated MW:

19kDa/26kDa/198kDa/209kDa/215kDa

Applications:

WB

Reactivity:

Mouse

Antibody Information

Recommended dilutions:

WB 1:500 - 1:1000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

4065

Uniprot

O60449

Synonyms:

LY75; CD205; CLEC13B; DEC-205; GP200-MR6; LY-75


Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 440-680 of human LY75 (NP_002340.2).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot analysis of extracts of mouse lung, using LY75 antibody (CAB10065) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (CABS014) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit (CABM00020). Exposure time: 90s.