

Product Information

Product SKU:	CAB10148	Gene ID:	55897	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	20kDa	Conjugate:	Unconjugated
Calculated MW:	29kDa	Isotype:	IgG

Immunogen Information

**Background:** Enables DNA-binding transcription factor activity and transcription cis-regulatory region binding activity. Involved in several processes, including endothelial cell differentiation; heart development; and positive regulation of transcription by RNA polymerase II. Predicted to be part of chromatin. Predicted to be active in nucleus.

Recommended Dilution:	WB,1:200 - 1:1000 IHC-P,1:500 - 1:1000
Synonyms:	bHLHc5; MESP1
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 149-268 of human MESP1 (NP_061140.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.