

CAB10236

Product Information

Product SKU:	CAB10236	Gene ID:	5331	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	160kDa	Conjugate:	Unconjugated
Calculated MW:	139kDa	Isotype:	IgG

Immunogen Information

Background: This gene encodes a member of the phosphoinositide phospholipase C beta enzyme family that catalyze the production of the secondary messengers diacylglycerol and inositol 1,4,5-triphosphate from phosphatidylinositol in G-protein-linked receptor-mediated signal transduction. Alternative splicing results in multiple transcript variants.

Recommended Dilution: WB,1:500 - 1:2000 IHC-P,1:50 - 1:200

Synonyms: SMDCD; PLC beta 3 (PLCB3)

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 1010-1234 of human PLC beta 3 (PLCB3) (NP_001303243.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.