

TUBGCP2 Rabbit Polyclonal Antibody



CAB10283

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

102kDa

Calculated MW:

88kDa/102kDa/105kDa

Applications:

WB

Reactivity:

Human, Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

10844

Uniprot

Q9BSJ2

Synonyms:

TUBGCP2; GCP-2; GCP2; Grip103; SPBC97; Spc97p; h103p; hGCP2; hSpc97


Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 1-180 of human TUBGCP2 (NP_006650.1).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot analysis of extracts of various cell lines, using TUBGCP2 antibody (CAB10283) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (CABS014) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit (CABM00020). Exposure time: 5s.