

CAB10321

Product Information

Product SKU:	CAB10321	Gene ID:	64101	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	73kDa	Conjugate:	Unconjugated
Calculated MW:	73kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to be involved in modulation of chemical synaptic transmission and synapse organization. Predicted to act upstream of or within synapse organization. Predicted to be located in neuron spine and postsynaptic membrane. Predicted to be active in Schaffer collateral - CA1 synapse; glutamatergic synapse; and plasma membrane. Predicted to be integral component of postsynaptic density membrane.

Recommended Dilution:	WB,1:1000 - 1:2000
Synonyms:	NAG14; NGL-2; LRRC4
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 544-653 of human LRRC4 (NP_071426.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.