

## CAB10423

---

### Product Information

<b>Product SKU:</b>	CAB10423	<b>Gene ID:</b>	133418	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	-	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human

---

### Additional Information

<b>Observed MW:</b>	37kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	37kDa	<b>Isotype:</b>	IgG

---

### Immunogen Information

**Background:** This gene encodes a transmembrane glycoprotein that is a member of the immunoglobulin superfamily. The encoded protein may be involved in cell growth and development by mediating interactions between the cell and extracellular matrix. A pseudogene of this gene is found on chromosome 1.

<b>Recommended Dilution:</b>	WB, 1:500 - 1:2000
<b>Synonyms:</b>	GP70; EMB
<b>Purification Method:</b>	Affinity purification
<b>Immunogen:</b>	Recombinant fusion protein containing a sequence corresponding to amino acids 33-260 of human EMB (NP_940851.1).
<b>Storage:</b>	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.