

# ATF7IP Rabbit Polyclonal Antibody



CAB10446

---

## Product Information

### Size:

20uL, 50uL, 100uL, 200uL

### Observed MW:

### Calculated MW:

118kDa/136kDa/137kDa

### Applications:

IHC

### Reactivity:

Human

## Antibody Information

### Recommended dilutions:

IHC 1:50 - 1:200

### Source:

Rabbit

### Isotype:

IgG

### Purification:

Affinity purification

## Protein Background

### Immunogen information

#### Gene ID:

55729

#### Uniprot

Q6VMQ6

#### Synonyms:

ATF7IP; AM; ATF-IP; MCAF; MCAF1; p621

#### Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 525-725 of human ATF7IP (NP\_060649.3).

#### Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---