

Product Information

Product SKU:	CAB10446	Gene ID:	55729	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse

Additional Information

Observed MW:	250kDa	Conjugate:	Unconjugated
Calculated MW:	136kDa	Isotype:	IgG

Immunogen Information

Background: ATF7IP is a multifunctional nuclear protein that associates with heterochromatin. It can act as a transcriptional coactivator or corepressor depending upon its binding partners (summary by Liu et al., 2009 [PubMed 19106100]).

Recommended Dilution:	WB,1:500 - 1:1000
Synonyms:	AM; MCAF; p621; MCAF1; ATF-IP; ATF7IP1; ATF7IP
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 525-725 of human ATF7IP (NP_060649.3).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.