

CAB10498

Product Information

Product SKU:	CAB10498	Gene ID:	55107	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse

Additional Information

Observed MW:	150kDa	Conjugate:	Unconjugated
Calculated MW:	114kDa	Isotype:	IgG

Immunogen Information

Background:	Enables calcium activated cation channel activity; intracellular calcium activated chloride channel activity; and iodide transmembrane transporter activity. Involved in cation transport; inorganic anion transport; and positive regulation of insulin secretion involved in cellular response to glucose stimulus. Located in apical plasma membrane and nucleoplasm.
Recommended Dilution:	WB,1:1000 - 1:2000 IF/ICC,1:100 - 1:200
Synonyms:	DOG1; INDMS; TAOS2; ORAOV2; TMEM16A; ANO1
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 400-500 of human ANO1 (NP_060513.5).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.